

iVascular[®]

therapies for living


luminor 14m

Paclitaxel eluting PTA balloon dilatation catheter

Safe and precise
delivery of Paclitaxel

LUMINOR Registry
207 patients

Nanotechnology
dosage system.
Ultra-thin coating



Assured coating
stability during
navigation



Proprietary coating
technology **TransferTech**



Homogeneous
& **anti-cracking** coating



www.ivascular.global

luminor 14m features

- > **Over-the-wire catheter (OTW)**
- > Usable catheter length: **100 cm** and **150 cm**
- > Recommended guidewire: **0.014"**
- > Tip profile: **0.017"**
- > Crossing profile: **from 0.026" up to 0.042"**
- > Nominal pressure: **7 atm**
- > RBP: **16 atm**
- > ABP: **22 atm**

Product with CE mark, certified by Notified Body 0318

Usable catheter length	Balloon Ø (mm)	Balloon Length (mm)					
		40	60	80	120	150	200
100 cm	1.5	BP DPC14M 100 150 040	BP DPC14M 100 150 060	BP DPC14M 100 150 080	-	-	-
	2.0	BP DPC14M 100 200 040	BP DPC14M 100 200 060	BP DPC14M 100 200 080	BP DPC14M 100 200 120	BP DPC14M 100 200 150	BP DPC14M 100 200 200
	2.5	BP DPC14M 100 250 040	BP DPC14M 100 250 060	BP DPC14M 100 250 080	BP DPC14M 100 250 120	BP DPC14M 100 250 150	BP DPC14M 100 250 200
	3.0	BP DPC14M 100 300 040	BP DPC14M 100 300 060	BP DPC14M 100 300 080	BP DPC14M 100 300 120	BP DPC14M 100 300 150	BP DPC14M 100 300 200
	3.5	BP DPC14M 100 350 040	BP DPC14M 100 350 060	BP DPC14M 100 350 080	BP DPC14M 100 350 120	BP DPC14M 100 350 150	BP DPC14M 100 350 200
	4.0	BP DPC14M 100 400 040	BP DPC14M 100 400 060	BP DPC14M 100 400 080	BP DPC14M 100 400 120	-	-
150 cm	1.5	BP DPC14M 150 150 040	BP DPC14M 150 150 060	BP DPC14M 150 150 080	-	-	-
	2.0	BP DPC14M 150 200 040	BP DPC14M 150 200 060	BP DPC14M 150 200 080	BP DPC14M 150 200 120	BP DPC14M 150 200 150	BP DPC14M 150 200 200
	2.5	BP DPC14M 150 250 040	BP DPC14M 150 250 060	BP DPC14M 150 250 080	BP DPC14M 150 250 120	BP DPC14M 150 250 150	BP DPC14M 150 250 200
	3.0	BP DPC14M 150 300 040	BP DPC14M 150 300 060	BP DPC14M 150 300 080	BP DPC14M 150 300 120	BP DPC14M 150 300 150	BP DPC14M 150 300 200
	3.5	BP DPC14M 150 350 040	BP DPC14M 150 350 060	BP DPC14M 150 350 080	BP DPC14M 150 350 120	BP DPC14M 150 350 150	BP DPC14M 150 350 200
	4.0	BP DPC14M 150 400 040	BP DPC14M 150 400 060	BP DPC14M 150 400 080	BP DPC14M 150 400 120	-	-

● Introducer compatibility: **4F**

Distributed by:



Manufactured by:

Life Vascular Devices Biotech, S.L.
www.ivascular.global
info@ivascular.global